
#### One vision, Two sharp eyes with Our Innovation

# AL-100

### **Biometer**

Capturing Critical Eye Measurements Maximum Precision and Ease-of-use







One vision, Two sharp eyes with Our Innovation

## AL-100

#### **Measure Axial Length Quickly & Easily**

The AL-100 is designed for fast, simple operation. Tone-assisted measurement notifies the operator when the probe is aligned & measurements recorded. The AL-100 takes up to 10 readings, for use in the IOL calculation. The AL-100 measures axial lengths from 15-40mm with an accuracy of  $\pm 0.1$ mm & a resolution of 0.01mm. Eye modes include Normal, Aphakic, Pseudophakic & Dense Cataract. Pseudophakic mode is able to select three materials:Silicone, Acrylic and PMMA. Measurements may be acquired automatically or manually & in Contact or Immersion modes.

#### **IOL Power Calculations**

The AL-100 calculates IOL power with a vast array of formulas: SRK II, SRK/T, Holladay, Haigis Standard, Haigis Optimized, Showa, Hoffer®-Q & Shammas-PL. Other formulas may be added as well. The AL-100 can display up to two lens constants & the corresponding IOL powers. Up to 10 lenses & lens constants can also be stored for use (along with the corresponding Surgeon Factor, SF).

#### **Data Storage & Output**

Entering patient & lens data is simple & direct with the large, easy-to-read touch screen. No more small or unattached keypads to make data entry a chore. The AL-100 also stores up to six patients in its internal memory & up to 240 patients in the optional 512Kb PC memory card. The AL-100 also has built-in printer, which can print either a standard waveform & IOL calculation, for both eyes, or more complete printout with all of the measured & calculated data, for both eyes.

In addition, there is an RS-232C port for direct connection to a PC & the office database.

#### **AL-100 SPECIFICATIONS**

#### **Biometry / IOL Power Calculation**

**Measurement Range** 

Axial Length 15.00 – 40.00mm
Anterior Chamber Depth 1.8 – 7.0mm
Lens Thickness 2.0 – 6.0mm

**Accuracy** 

Measurement Accuracy ± 0.1mm Resolution ± 0.01mm

**IOL Calculation** 

Formula SRK II, SRK/T, HOLLADAY, SHOWA, HAIGIS Optimized, HAIGIS Standard

Hoffer®-Q, Shammas-PL


**Biometry Probe** 

Type Solid (with Immersion option)

Fixation Light Built-in (Red LED)

Ultrasound Frequency 10 MHz Diameter of the Tip 5.3mm

Size / Weight 8 mm x 97 mm / 30 g (0.3 x 3.8'' / 1 oz.)



Clear spikes for accurate & reliable measurements.

Phy:			Utility	14 02
No	AXIAL	ACD	LENS	CSI .
1	25.64	3.78	4.18	C D H
2	25.64	3.78	4.17	SelAvg
3	25.64	3.78	4.18	Selling
4	25.65	3.81	4.16	A -
Avg★	25.64	3.79	4.17	

Quick & easy date review.



Easy comparison with different parameters.

5.7 inches Color 220mm (W) x 222mm (D) x 275mm (H) 8.7" (W) x 8.7" (D) x 10.8" (H) 4kg (8.8lbs.)

AC 100 to 240V 50/60 Hz 35/42VA

### Weight Power Source Voltage Frequency

**Main Body** 

**Dimensions** 

Power Consumption

STN Liquid Crystal Display



#### **Tomey Corporation [Asia-Pacific]**

2-11-33 Noritakeshinmachi Nishi-ku, Nagoya, 451-0051, Japan Tel: ++81-52-581-5327 Fax: ++81-52-561-4735 E-mail: intl@tomey.co.jp

For more information, visit our web site

#### Tomey GmbH [Europe]

Wiesbadener Straße 21 90427 Nürnberg, Germany Tel: ++49-911-9385462-0 Fax: ++49-911-9385462-20 E-mail: info@tomey.de

http://www.tomey.com